

FORM 11 (ENG): ACCREDITED CALIBRATION - Order & Decontamination Declaration

Customer/no.: _____ Contact (name): _____
 Address: _____ Department: _____
 _____ E-mail: _____
 EAN no.: _____ Phone: _____
 Service agreement no.: _____ PO/order: _____

No.	Pipette manufacturer:	Serial no.:	Customer ID no.:	No.	Pipette manufacturer:	Serial no.:	Customer ID no.:
1				6			
2				7			
3				8			
4				9			
5				10			

	5.3	6.3	10.3	10.1	10.1 Fixed vol.	Pipetting mode - write P, rP or D If the box is empty; Pipettes will be calibrated in P-mode and Dispensers in D-mode	
PreCal (As found - before service)*						**Multichannel pipettes If no X - all channels will be calibrated	
Calibration (As left - after service)						All channels are be calibrated	
ONLY calibration (No service)						3 channels are calibrated The outer channels and the middle one	

*) A pre-calibration documents the functionality of the pipette before a service is done and is recommended according to ISO 17025 requirements.
 **) For multichannel pipettes where one or more calibration points are $\leq 10 \mu\text{l}$, Dandiag uses a 5-decimal multichannel balance instead of a 6-decimal balance as otherwise prescribed in ISO 8655-6. When using a multichannel balance, use of the multichannel pipette is reflected and the uncertainty is approx. $0,07 \mu\text{l}$ against approx. $0,04 \mu\text{l}$ for a 6-decimal balance.

SPECIFICATIONS	TIPS / Tip cone FILTERS	SERVICE, ADJUSTMENT etc.
If no X - the ISO 8655 max. permissible errors will be used		
ISO 8655 - Max. errors	Customers tips - should be sent in original packaging	Service is ordered
TT Type test - manufac. performance spec.	NO thanks to inserting filter	NO adjustment
TT*1,5 Can be use for not new pipettes		Repair, please write the malfunction
Customers own user specifications must be attached		Pipettes are purchased with an annual service concept

Remarks: _____

Please contact me, if the repair price per pipette exceeds: _____ DKK (don't include change of battery)

The undersigned hereby declares that the pipettes have been carefully cleaned, decontaminated and **filters are removed**:

70% Ethanol Decontamination liquid Other

The pipettes have been use for?

The pipettes have been used for radioactive liquids

If the pipette has been used for radioactive liquids, documentation for residual radiation must be sent (attached).

Date: _____ Name: _____